

Flameproof (Spark-free) Walk-In Chambers



The name you can trust

FIREsafe



For Safety & Security of your volatile, flammable materials

Specifically designed to store aerosols, volatile, flammable liquids, reagents, chemicals and organic solvents.



Fire proof Panels



Ethernet Connectivity



21 CFR Part 11
Compliance Software



Standby Cooling Systems



Mobile Alert

FIREsafe

Introducing Flameproof (Spark-free) Walk-in Chambers

Features

- ◆ FM approved PUF Panels
- ◆ Spark-free Interior & Exterior
- ◆ Flame-proof Components
 - Motors | HMI | Light Fittings | Control Panel | Sensors*
- ◆ FM200 Fire Suppression System
 - Automatic Start Function | AV Alarms/Hooters | Flame & Smoke Detection Sensors
- ◆ Complete with Validation & Test Certificates
- ◆ Safety features with comprehensive alarms



Safety Features: FM200 Fire Suppression System

- ◆ In case of fire inside the chamber, Inert gas will be purged
- ◆ Automatic Start Function
- ◆ Audio-Visual Alarms/Hooters
- ◆ Flame & Smoke Detection Sensors installed inside the Walk-In Chamber
- ◆ Supply to all instrumentation will be cut-off
- ◆ No risk of personal injury
- ◆ Faster than other extinguishing systems using natural gases
- ◆ Extinguishing agent storage takes up little space as consumption is very low
- ◆ Electrically non – conducting
- ◆ Best environmental compatibility of the HFCs
- ◆ No clean up required after discharge
- ◆ Modular design using multiplex extinguishing agent storage containers
- ◆ Other Walk-In Chambers also available in Flameproof



Our Products Range: Humidity (Stability) Chambers | Photostability Chambers | Freeze Thaw Chambers | Cooling Chambers | Walk-in Chambers | Deep Freezers | Plasma Freezers(-40°C) | Ultra Low Temperature Freezers(-80°C) | Incubators

ISO 9001:2008 Certified Company



The name you can trust

Visit us at : www.newtronic.in

NEWTRONIC LIFECARE EQUIPMENT PVT. LTD.

H.O.:108-ABCD, Kandivali Co-Op. Ind. Estate Ltd.
Charkop, Kandivali (West), Mumbai- 400067, India
Tel: +91-22-2867 9326 | Fax: +91-22-2867 6059

Plant: 1205-1206, GIDC, Road No.13
At - Umbergaon, Dist. Valsad - 396171, Gujarat

Email:

India : sales@newtronic.in
International : exports@newtronic.in

Service:

service@newtronic.in
service.exports@newtronic.in