



The name you can trust

Blood Bank Refrigerator

Temp. Range: 2°C to 6°C



Leading Solution For Blood Bank Refrigeration

The blood bank refrigerator is an essential piece of equipment in the immuno- hematology department and provides safe and convenient storage of whole blood, blood components (e.g, blood cells, plasma), and reagents. Blood bank refrigerators ensure freshness and integrity of blood and blood components. Our Blood Bank Refrigerators meet the highest regulatory requirements.

Ergonomically designed chambers
for single hand operations



Cyclic Defrost



Forced Air Draft
Circulation



Uniform
Temperature



Low Noise



Environment
Friendly



Intelligent Alarms

Blood Bank Refrigerator

Temp. Range: 2°C to 6°C

Standard Features

Inner body	Ss304
Outer body	Pre Coated Galvanised Iron
Insulation	High density CFC free PUF insulation
Door	Multi layer transparent door with lock & handle Magnetic gasket is provided for airtight sealing
Moveability	Castor wheels are provided for easy movement with an option for having locks
Controller	Electronic controller with digital display
Lighting	Bright LED interior lighting
Alarms	Temperature deviation alarm and door ajar alarm
Cooling System	CFC free cooling system. The compressor is hermetically sealed for better efficiency It has low noise and has low power consumption.
Power supply	230 Volts 50 Hz main supply. Also available with 60 Hz cycle
Safety	Anti-frost safety device.



Selection Guide

Model	Capacity		Internal Size (H*W*D) mm	External Size (H*W*D) mm	Drawers Qty.	Single Drawer Blood Bag Capacity	All Drawer Blood Bag Capacity	Equipment Weight (kg)
	CuFt.	Liter						
NLBR7MI	7	200	825 X 490 X 495	1445 X 590 X 610	3	36	108	150
NLBR11MI	11	300	1000 X 600 X 500	1675 X 700 X 615	4	52	208	200
NLBR14MI	14	400	1000 X 600 X 670	1675 X 700 X 785	4	76	304	250
NLBR17MI	17	500	1200 X 600 X 695	1875 X 700 X 810	5	76	380	300
NLBR21MI	21	600	1255 X 650 X 775	1945 X 750 X 890	5	76	380	350


Optional Accessories

Data Logger To automatically record parameters

Circular chart recorder To archive data points over a time interval



Data Logger




Circular Chart Recorder




Adjustable Trays



Digital Controller with Audio Visual Alarm



Bright LED Lights



Low Power Consumption



Low Noise Cooling System



India : sales@newtronic.in
International : exports@newtronic.in
Visit us on: www.newtronic.in

NEWTRONIC LIFECARE EQUIPMENT PVT. LTD.
ISO 9001:2008 Certified Company Safety Gloves

H.O.: Kanakia Wall Street, A-Wing 612, 6th Floor, Chakala, Andheri-Kurla Road, Andheri (E), Mumbai-400 093, India. T: +91-22-6853 4600

Plant: Plot No-1205 & 1206, G.I.D.C. Ind. Estate, Road No.: 13, At Umbergaon – 396171, Dist.-Valsad (Gujarat). India.

